

SKINCOOL ICE ROLLER


Freeze it Roll it Enjoy it

Product Description:

Ice roller is a hand held device that has a water and gel filled rolling head that can be effectively frozen in the freezer. It can be used at any place and any time to release pain and discomfort instantly. Ice roller can be also used daily to revitalize, refresh and rejuvenate the skin. It is a simple and easy device that makes your life convenient without melting ice.

Functions:

1. Ice Roller can be also used daily to revitalize, refresh and rejuvenate the skin. It is a simple and easy device that makes your life convenient without melting ice.
2. Ice Roller can be use before and after IPL/laser, botox, micro needle therapy and microdermabrasion treatment and any cosmetic procedure. It minimizes discomforts, irritation and redness.
3. Use Ice roller on your face, neck and chest every morning (Ice cube facials). Your face will look revitalized and young. It also prevents wrinkle formation and lifts your face.
4. Use Ice roller on your eyes when you feel fatigue. It refreshes. Roll over on any parts of your body for instant pain relief. It is a must have item at home for sport injury, burn and insect bites without making mess with melted ice. After shaving, it will calm and soothe your skin.
5. Keep Ice Roller on fever for cooling effect.



How to use:

1. Clean the roller head and put in the bag
2. Keep Ice Roller in the freezer at all times or have a quick-freezing with 10-15 minutes.
3. Take the roller head from freezer and install with the handle
4. Enjoy the cooling and calming massaging effects on your face and body

Features:

1. Easy to use, quick cooling, long duration, convenient to carry
2. Wide application, suitable for all ages.
3. Convenient, safe and health

Caution:

1. Do not use on wounded skin tissue
2. Store in a cool place
3. Do not heat this product